

GE Healthcare

# Dash 2500

The standard of excellence  
for lower-acuity monitoring





# The clinical

With monitored patients, even small measurement variations can indicate significant changes in your patient's condition. You want assurance that any deviation is related to the patient – not the monitor.

The Dash®2500 monitor from GE Healthcare allows you to deliver a new standard of clinical excellence to patients in lower-acuity settings such as ER, recovery/PACU, outpatient surgery, GI/endoscopy, progressive care, interventional radiology and special care units. As a member of the Dash family of monitors, the Dash 2500 leverages the same parameters and algorithms as GE's clinically advanced, higher-acuity monitors. With a consistent foundation of reliable and accurate clinical intelligence, you can confidently make better, faster decisions.

The Dash 2500 treats every patient like you do, with unwavering quality, speed and precision.

# standard of excellence

## **Rugged and appealing**

The Dash 2500 monitor is designed for demanding clinical environments. With outstanding durability, it meets the same rigorous requirements as all Dash monitors. Its chemically resistant casing is made of GE LEXAN® plastic, recently reformulated for superior resilience and enduring performance. So even tough conditions and powerful cleansers won't damage it.

Though built to last, the Dash 2500 is also ergonomic and easy to use. The display is easily configurable to suit your viewing preferences and information requirements. And its compact size and long battery life make it the perfect portable monitor; you can transport patients anywhere in the hospital – and maintain an uninterrupted flow of data.

With the Dash 2500 monitor, there's no need to sacrifice performance in the interest of cost. At GE Healthcare, we apply our most advanced technologies, expertise and legacy of innovation to every monitor. It's part of our vision to unite the way people, devices and technology work together to create a fluent system of clinical intelligence. So you always have the information you need to deliver the best, most cost-efficient care – everywhere.

# technology of today

The Dash 2500 treats every patient like you do, with unwavering quality, speed and precision.

## **Robust and economical**

The Dash 2500 is a full-featured, affordable bedside monitor. Its powerful clinical parameters capture vital patient measurements while unparalleled arrhythmia detection algorithms help to accurately detect lethal arrhythmias. With Dash 2500 by your side, you can spend less time chasing false alarms and more time responding to your patients' true care needs.

The Dash 2500 includes:

- EK-Pro™ arrhythmia program
- DINAMAP® SuperSTAT™ non-invasive blood pressure
- Masimo SET® or Nellcor® OxiMax® SpO2
- Alaris Turbo Temp®

The Dash 2500 gives you the clinical intelligence you need to manage your patients vigilantly, intervene appropriately and return them to health as soon as possible.





©2007 General Electric Company.

GE and GE Monogram are trademarks of General Electric Company.  
Masimo and SET are trademarks of Masimo Corporation.

Nellcor and OxiMax are trademarks of Nellcor Puritan Bennett, Inc.

Turbo Temp is a trademark of Cardinal Health, Inc.

Dash, DINAMAP, EK-Pro, LEXAN and SuperSTAT are trademarks of GE Medical Systems, Information Technologies, Inc.

General Electric Company, doing business as GE Healthcare.

## Healthcare Re-imagined

GE is dedicated to helping you transform healthcare delivery by driving critical breakthroughs in biology and technology. Our expertise in medical imaging and information technologies, medical diagnostics, patient monitoring systems, drug discovery, and biopharmaceutical manufacturing technologies is enabling healthcare professionals around the world discover new ways to predict, diagnose and treat disease earlier. We call this model of care "Early Health." The goal: to help clinicians detect disease earlier, access more information and intervene earlier with more targeted treatments, so they can help their patients live their lives to the fullest. Re-think, Re-discover, Re-invent, Re-imagine.

GE Healthcare

P.O. Box 900, FIN-00031 GE, Finland

Tel. +358 10 394 11 • Fax +358 9 146 3310

[www.gehealthcare.com](http://www.gehealthcare.com)



GE imagination at work