



M · O · W
Master Of World Medical

MOW Silicone Resuscitator

Made of special silicone and medical grade plastic

- The valve made of special silicone had no resistance.
- You can use a autoclave all plastic parts.

Many disinfects are possible

- Sterilization in low temperature.
- Disinfect by bolling / gas / freezing / steam up to 136'c/ hight voltage.

New type mask

- This new type mask doesn't need air possible.

Economical Efficiency.

- Both of the replacement and the repair for broken parts are possible.

Oxygen supply is 100% possible

- Oxygen supply is 100% possible by attaching reservoir valves and bag.
- And the most effects are generated by the least oxygen.

Easy to regardless and operation

- The pressure bag of silicone can be re-expanded easily regardless of extreme temperature,
- And the special treatment of surface makes rhe operation so easy

MOW Silicone Resuscitator



(Adult) : weight 30Kg more than



(Child) : weight 7~30Kg



(Infant) : weight 7Kg and less

Oxygen Reservoir



(Adult/Child) : 2.6ℓ



(Infant) : 0.6ℓ



Air Way
(No.40mm – 110mm.)

Resuscitation Mask





M · O · W
Master Of World Medical

ชุดช่วยหายใจแบบมือบีบ

ลูกยางบีบช่วยหายใจ (Silicone Bag) ทำจากซิลิโคนแบบพิเศษและพลาสติกอย่างดีที่ใช้ในการแพทย์ อุปกรณ์ทุกส่วนสามารถถอดล้างทำความสะอาด ด้วยวิธี Autoclave ได้



* (ผู้ใหญ่) : น้ำหนัก 30 กิโลกรัมขึ้นไป



* (เด็กโต) : น้ำหนัก 7~30 กิโลกรัม



* (เด็กเล็ก) : น้ำหนักต่ำกว่า 7 กิโลกรัม

มีมือจับเพื่อความสะดวกในการใช้



ถุงสำรองออกซิเจน (Oxygen Reservoir) สามารถสำรองอากาศได้ 100%



(ผู้ใหญ่และเด็กโต) : 2.6ลิตร



(เด็กเล็ก) : 0.6ลิตร

หน้ากาก (Resuscitation Mask) ทำจากซิลิโคนใสทั้งอัน



Air Way



(No.40mm – 110mm.)