

Description:



A-100B-S is a Time Fetal Doppler, which is designed for clinic and small hospital. It is a lightweight and space-saving Time Fetal Doppler, which can be operated at bedside and acted as handset. It is designed for use in the antepartum applications.

Features:

- Indicates FHR as numeral and outputs fetal heart sound.
- Patented acoustic design enhanced vivid sound.
- High Sensitive 2.5 MHz Doppler probe.
- RTC Timer inside can display day/hour/minute.
- Built-in rechargeable battery for 24 hours continuous work
- Medical-use switch power supply guarantee 50,000hrs working life.

Specifications:

FHR:

- Ultrasound Frequency: 2.5MHZ
- Intensity : 2.5 mW/cm² or less
- FHR Range:60-210BPM
- FHR Accuracy: ±1BPM

Sound:

- Doppler Sound with volume control
- Alarm Sound for improper FHR

Display:

- 6-Segment LED to display FHR or Hour/Minute /Second
- Battery capacity display
- Power off/on display
- Alarm on/off display
- Probe off display

Function:

- Alarm on/off function
- Volume +/- function
- Menu set function